

Product Name :
Orthopedic Bed, Fowler**Product Code :**
NMCBS106004**Description :**

Orthopedic Bed, Fowler
Orthopaedic bed with mechanically operated back rest 0 to 70° & knee rest 0 to 35° all by Super Smooth crank mechanism. Size: 2150 x 900 x 500 mm. Epoxy Coated Steel frame work and 4 section perforated top. Stainless Steel Orthopaedic attachment (without weight chain and weight hangers). Head panel and foot panel equipped with safety lock and roller bumpers. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction
Orthopedic Bed, Fowler,
Hospital Bed exporters, Hospital Bed manufacturers