

**Product Name :**  
Birthing Bed Motorized**Product Code :**  
NMCBS106012**Description :**

Birthing Bed Motorized 1. Overall approximate dimensions 2050 mm (l) x 1020 mm (w) 2. Minimum height 490 mm (without mattress) 3. Maximum height 780 mm (without mattress) 4. Detachable abs head and foot boards. 5. Collapsible abs guardrails. 6. 125mm double wheel castors (two with brakes) 7. Electrical lifting system with (taiwan) motors. (5 year warranty) 8. Bed framework made from mild steel erw rectangular tubes. 9. Bed top made from nonmagnetic bakelite sheet. X-ray compatible. 10. Trendelenburg & reverse trendelenburg. 11. Telescopic leg section that reduces the overall length of the bed, allowing maximum access to mother & baby. 12. Comprehensive range of easy to fit accessories including stirrups leg rests, hand grips, high and low foot supports for accessories step and a fluid collection bowl. 13. Bright, modern appearance. Birthing Bed Motorized, Hospital Bed exporters, Hospital Bed manufacturers