

Product Name :
Naugramedical Fowler Bed Manual**Product Code :**
NMCBS106013**Description :****Naugramedical Fowler Bed Manual India China****Technical Specification :**

Naugramedical Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism.

Size: 2140 940 500 mm.

Two easy lifting guardrails (1 on each side), which are safe and reliable, and can be fixed upward and downward.

Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel.

Epoxy coated mild steel frame work and 4 section perforated top.

125mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press.

Provision for I.V. Rod. Fowler Bed Manual, Hospital Bed exporters, Hospital Bed manufacturers

Medical Hospital Equipments Manufacturers

Website: [SITE](#).

NAUGRA[®]