

Product Name :
Fowler Steel Bed Manual 2 Function**Product Code :**
NMCBS106015**Description :**

Fowler Steel Bed Manual 2 Function

Technical Specification :

Fowler Steel Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism.

Size: 2150 98 500 mm.

Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel.

Collapsible Aluminium side railings.

Equipped with Food Tray & with shelf at the bottom for keeping blankets.

Epoxy coated mild steel frame work and 4 section perforated top.

125mm dia noiseless castors (all with individual brakes).

Provision for I.V. Rod on both sides of the bed.

Freight saving knock-down construction.

Fowler Steel Bed Manual 2 Function, Hospital Bed exporters, Hospital Bed manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,