

Product Name :
Orthopaedic Bed Two Sectioned Adjustable Knee Rest**Product Code :**
NMCBS106017**Description :**

Orthopaedic Bed Two Sectioned Adjustable Knee Rest
• Orthopaedic bed with mechanically operated back rest & knee rest all by crank mechanism.
• Knee-rest portion sectioned in two parts for patient's convenience.
• Size: 2030 x 900 x 500 mm.
• Epoxy coated mild steel 5 section top.
• Frame work made of stainless steel.
• Balkan Frame with monkey bars and monkey chains.
• Pulley attachment (without weight chain and weight hangers).
• Mounted on protective stumps.
• Provision for I.V. Rod.
• Freight saving knock-down construction.
Type:- Without castor With castor Orthopaedic Bed Two Sectioned Adjustable Knee Rest, Hospital Bed exporters, Hospital Bed manufacturers