

**Product Name :**  
Medical Semi Fowler Beds**Product Code :**  
NMCBS106018**Description :****Hospital Semi Fowler Beds, Medical Semi Fowler Beds India China****Technical Specification :**

Semi Fowler Bed Is Backrest section maneuvered by crank mechanism to obtained requisite position.

Semi Fowler Bed is Framework of precise rectangular steel tube mounted on PVC stumps.

Semi Fowler Bed Has two sections of perforated CRCA sheets duly machine press bent.

Semi Fowler Bed Tubular head and foot bows unequal heights with vertical tubes supports on head bow.

Semi Fowler Bed -Provisions for four saline rod locations

Semi Fowler Bed is Overall Size: 198L x 90W x 60 H cms

Finish: Pre-treated to prevent rust and Epoxy Powder Coated. Hospital Semi Fowler Beds, Hospital Bed exporters, Hospital Bed manufacturers

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**