

Product Name :
Electrical Intensive Care Bed**Product Code :**
NMCBS106022**Description :**

Electrical Intensive Care Bed Intensive care Hospital Bed for comfortable recovery of patient under cumbersome conditions. $\dot{\imath}\dot{\imath}\frac{1}{2}$ Four section CRCA sheet Top, uniformly perforated fixed to the frame work & base made of ERW Steel tube /rectangular tubes $\dot{\imath}\dot{\imath}\frac{1}{2}$ MS channels mounted on four special Swivelling castors with 2 brakes $\dot{\imath}\dot{\imath}\frac{1}{2}$ full length collapsible stainless steel safety railing. $\dot{\imath}\dot{\imath}\frac{1}{2}$ Adjustable telescopic I.V.Rod. Remote control system offers:

- $\dot{\imath}\dot{\imath}\frac{1}{2}$: Electrical Intensive Care Bed Is High & low position
- $\dot{\imath}\dot{\imath}\frac{1}{2}$: Electrical Intensive Care Bed is Headrest position
- $\dot{\imath}\dot{\imath}\frac{1}{2}$: Electrical Intensive Care Bed Trend lenburg/reverse trend lenburg position
- $\dot{\imath}\dot{\imath}\frac{1}{2}$: Electrical Intensive Care Bed Size: 198(L) x 90(W) x 60 to 80(H) cms
- $\dot{\imath}\dot{\imath}\frac{1}{2}$ Finish: Pre-treated to prevent rust and Epoxy Powder Coated.

Electrical Intensive Care Bed, Hospital Bed exporters, Hospital Bed manufacturers