

**Product Name :**  
Electrical Intensive Care Bed**Product Code :**  
NMCBS106022**Description :**

Electrical Intensive Care Bed Intensive care Hospital Bed for comfortable recovery of patient under cumbersome conditions.  $\dot{\imath}\dot{\imath}\frac{1}{2}$  Four section CRCA sheet Top, uniformly perforated fixed to the frame work & base made of ERW Steel tube /rectangular tubes  $\dot{\imath}\dot{\imath}\frac{1}{2}$  MS channels mounted on four special Swivelling castors with 2 brakes  $\dot{\imath}\dot{\imath}\frac{1}{2}$  full length collapsible stainless steel safety railing.  $\dot{\imath}\dot{\imath}\frac{1}{2}$  Adjustable telescopic I.V.Rod. Remote control system offers:

- $\dot{\imath}\dot{\imath}\frac{1}{2}$  : Electrical Intensive Care Bed Is High & low position
- $\dot{\imath}\dot{\imath}\frac{1}{2}$  : Electrical Intensive Care Bed is Headrest position
- $\dot{\imath}\dot{\imath}\frac{1}{2}$  : Electrical Intensive Care Bed Trend lenburg/reverse trend lenburg position
- $\dot{\imath}\dot{\imath}\frac{1}{2}$  : Electrical Intensive Care Bed Size: 198(L) x 90(W) x 60 to 80(H) cms
- $\dot{\imath}\dot{\imath}\frac{1}{2}$  Finish: Pre-treated to prevent rust and Epoxy Powder Coated.

Electrical Intensive Care Bed, Hospital Bed exporters, Hospital Bed manufacturers