

**Product Name :**  
Electrical Intensive Care Bed (With Remote)**Product Code :**  
NMCBS106023**Description :**

Electrical Intensive Care Bed (With Remote) : Electrical Intensive Care Bed (With remote) Intensive care Hospital Bed for comfortable recovery of patient under cumbersome conditions. : Electrical Intensive Care Bed (With remote) Four section CRCA sheet Top,uniformly perforated fixed to the frame work & base made of ERW Steel tube /rectangular tubes : Electrical Intensive Care Bed (With remote) MS channels mounted on four special Swivelling castors with 2 brakes : Electrical Intensive Care Bed (With remote) full length collapsible stainless steel safety railing. : Electrical Intensive Care Bed (With remote) Adjustable telescopic I.V.Rod. Remote control system offers: : Electrical Intensive Care Bed (With remote) High & low position : Electrical Intensive Care Bed (With remote) Headrest position : Electrical Intensive Care Bed (With remote) Trendlenburg/reverse trendlenburg position : Electrical Intensive Care Bed (With remote) Size: 198(L) x 90(W) x 60 to 80(H) cms : Electrical Intensive Care Bed (With Remote), Hospital Bed exporters, Hospital Bed manufacturers