

Product Name :
Electrical Intensive Care Bed (Mechanically Operated)**Product Code :**
NMCBS106024**Description :**

Electrical Intensive Care Bed (Mechanically Operated) is an intensive care hospital bed for comfortable recovery of patient under cumbersome conditions. This intensive care bed is available in various important positions as per doctor's regiments. The intensive care bed base & main frame work is of M.S Rectangular tube / MS channels . Four section M.S. Top uniformly perforated intensive care bed is fixed to the main frame . Backrest, knee rest & trendlenburg / reverse trendlenburg positions adjustable by separate screws lever mechanism. Swiveling castors, with two brakes. This intensive care bed is convenient transfer of patient from bed to trolley and vice-versa. This intensive care bed is SS detachable collapsible side railings provided. This intensive care bed is provision for four saline rod locations. This intensive care bed size: 198 (L) x 91 (W) x 60 (H) cms Finish : Pre-treated to prevent rust and Epoxy Powder Coated. Remote control system offers: This intensive care bed is High & low position This intensive care bed is Headrest position This intensive care bed Trendlenburg/reverse trendlenburg position This intensive care bed size: 198(L) x 90(W) x 60 to 80(H) cms Finish: Pre-treated to prevent rust and Epoxy Powder Coated Electrical Intensive Care Bed (Mechanically Operated), Hospital Bed exporters, Hospital Bed manufacturers