

**Product Name :**  
Electrical Intensive Care Bed (Mechanically Operated)**Product Code :**  
NMCBS106024**Description :**

Electrical Intensive Care Bed (Mechanically Operated) is a 1/2 Intensive care hospital bed for comfortable recovery of patient under cumbersome conditions. 1/2 Intensive Care Bed is Various important positions available as per Doctor's regiments. 1/2 Intensive Care Bed base & main frame work is of M.S Rectangular tube / MS channels . Four section M.S. Top uniformly perforated 1/2 Intensive Care Bed is fixed to the main frame . Backrest, knee rest & trendlenburg / reverse trendlenburg positions adjustable by separate screws lever mechanism. Swiveling castors, with two brakes. 1/2 Intensive Care Bed is Convenient transfer of patient from bed to trolley and vice-versa. 1/2 Intensive Care Bed is SS detachable collapsible side railings provided. 1/2 Intensive Care Bed is Intensive Care Bed is Provision for four saline rod locations. 1/2 Intensive Care Bed Size: 198 (L) x 91 (W) x 60 (H) cms 1/2 Finish : Pre-treated to prevent rust and Epoxy Powder Coated. Remote control system offers: 1/2 Intensive Care Bed is High & low position 1/2 Intensive Care Bed is Headrest position 1/2 Intensive Care Bed Trendlenburg/reverse trendlenburg position 1/2 Intensive Care Bed Size: 198(L) x 90(W) x 60 to 80(H) cms 1/2 Finish: Pre-treated to prevent rust and Epoxy Powder Coated Electrical Intensive Care Bed (Mechanically Operated), Hospital Bed exporters, Hospital Bed manufacturers