

Product Name :
Stretcher On Trolley**Product Code :**
NMCBS109006**Description :**

Stretcher On Trolley We are one of the leading organizations, engaged in offering our customers with Stretcher On Trolley. These are manufactured using premium raw material and modern machinery that are procured from the reliable vendors of the market. In addition to this, we have a team of expert quality controllers, which checks all the products on various parameters to assure flawlessness. Other Information: Overall Size: L 1830 mm x W 530 mm x H 870 mm Frame: 50 mm x 25 mm and CRCA tube Detachable Stainless steel top fixed on CRCA tube frame mounted on 150 mm dia. heavy duty wheels Provision for fixing I.V. rod Pre-treated and epoxy powder coated Stretcher On Trolley, Hospital Trolley exporters, Hospital Trolley manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,