



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Monitor Trolley For Scopic Surgery

**Product Code :**  
NMCBS109013



**Description :**

Monitor Trolley For Scopic Surgery Overall size W 660 mm x D 560 mm x H 1498 mm Features: Available in full Stainless steel & epoxy powder coated Adjustable Shelves 125 mm dia. heavy duty casters with brakes Monitor Trolley For Scopic Surgery, Hospital Trolley exporters, Hospital Trolley manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)