

**Product Name :**  
Syringe Pump Stand

**Product Code :**  
NMCBS109019



**Description :**

Syringe Pump Stand \* Syringe pump stand is with 4 adjustable shelves. \* These shelves helps to hold syringe pumps. \* With IV stand on top. \* Five 50mm double wheel castors with brake. \* Electrical connection is provided at the bottom of the stand. \* The Complete stand is mobile on 50mm dia twin castors Syringe Pump Stand, Hospital Trolley exporters, Hospital Trolley manufacturers