

Product Name :
Stretcher Trolley(S.S)**Product Code :**
NMCBS109021**Description :**

Stretcher Trolley(S.S) Stretcher Trolley(S.S) * Strong construction made from stainless steel round tubes. * Mounted on four metallic body castors 200mm wheel dia, two with brake. * Removable stretcher top made of stainless steel sheet supported on stainless steel tubular frame. * Supplied with Utility tray & oxygen cylinder cage. Stretcher Trolley(S.S), Hospital Trolley exporters, Hospital Trolley manufacturers