

**Product Name :**  
Emergency & Recovery Trolleys**Product Code :**  
NMCBS109027**Description :**

Emergency & Recovery Trolleys  
Emergency & Recovery Trolley Made of heavy duty steel tubes.  
Removable X-Ray translucent provided with a tray for holding X-ray cassettes underneath the entire length of the stretcher.  
Emergency & Recovery Trolley Height adjustment by hydraulic system, trendelenburg and reverse trendelenburg by pressure cylinder.  
Emergency & Recovery Trolley Fitted with 15cms dia castors-two with brakes and rubber buffers on all sides.  
Emergency & Recovery Trolley Provided with SS swing away railing, SS telescopic saline rod, storage tray and cylinder cage.  
Emergency & Recovery Trolley Size: Size: 210L x 65W x 65 to 95 H cms.  
Finish: Pre-treated to prevent rust & Epoxy powder coated.  
Emergency & Recovery Trolleys, Hospital Trolley exporters, Hospital Trolley manufacturers