

**Product Name :**  
Multipurpose Bed Side Lockers**Product Code :**  
NMCBS110005**Description :**

Multipurpose Bed Side Lockers We are engaged in offering our clients with Multipurpose Bed Side Lockers that are available in different sizes and designs at market leading prices. All the lockers are designed in compliance with international quality norms & standards, using premium raw material and latest technology. In addition to this, we have the capability to customize all the products as per the specifications detailed down by the customers. Other information:  $\frac{1}{2}$  Overall Size: W 410 mm X D 370 mm X H 900 mm  $\frac{1}{2}$  Size after packing: W 840 mm X D 500 mm X H 190 mm  $\frac{1}{2}$  Available in CKD form to reduce freight cost & local transportation  $\frac{1}{2}$  Complete ABS (Acrylonitrile Butadiene Styrene) plastic  $\frac{1}{2}$  Different colors to select only for bulk quantity)  $\frac{1}{2}$  OEM supplies available by incorporating your logo (only for bulk quantity)  $\frac{1}{2}$  Easy to assemble, no special tools required  $\frac{1}{2}$  Available with casters and without casters  $\frac{1}{2}$  Can be used in Hospitals / Hotels / Offices / Laboratory Multipurpose Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers