

Product Name :
BED SIDE LOCKERS**Product Code :**
NMCBS110007**Description :**

BED SIDE LOCKERS Bed side lockers made of high quality sheet metal. 1½ SS top. 1½ Hospital bed side lockers Provided with one storage cabinet. 1½ Bed side lockers optionally available with single draw under the top and space for keeping utilities. 1½ Standard dimension:160mm L x 160mm W x 320 mm H. BED SIDE LOCKERS, Hospital Lockers exporters, Hospital Lockers manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,