

Product Name :
ABS Bed Side Lockers**Product Code :**
NMCBS110008**Description :**

ABS Bed Side Lockers The Hospital ABS Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Bedside Locker has a stainless steel top and mounted on two swivel castors. ABS Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,