

**Product Name :**  
ABS Bed Side Lockers**Product Code :**  
NMCBS110008**Description :**

ABS Bed Side Lockers The Hospital ABS Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Bedside Locker has a stainless steel top and mounted on two swivel castors. ABS Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#)