

**Product Name :**  
Stainless Steel Surface Bed Side Lockers**Product Code :**  
NMCBS110010**Description :**

Stainless Steel Surface Bed Side Lockers The Hospital Stainless Steel Surface Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Stainless Steel Surface Bedside Locker has a stainless steel top and mounted on two swivel castors. Stainless Steel Surface Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,