

Product Name :
Infant Bed**Product Code :**
NMCBS114002**Description :**

Infant Bed We have established ourselves as a renowned manufacturer and supplier of Infant Bed. These products are widely appreciated by the clients for their longer functional life, low maintenance, easy handling and zero defects. All the products are designed and developed using premium raw material, which is sourced from the reliable vendors of the market. Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg & Reverse Trendelenburg 12 $\frac{1}{2}$): $\frac{1}{2}$ Infant Bed / Child Cot / Baby Bassinet, Size: 870mm L $\frac{1}{2}$ 530mm W $\frac{1}{2}$ 920-1070mm H $\frac{1}{2}$ Height Adjustable by Gas Spring 920-1070mm $\frac{1}{2}$ Molded Plastic Crib (without joints) with clear visibility for observer $\frac{1}{2}$ Smooth Rotating Casters 100mm Dia will provide silent transportation $\frac{1}{2}$ Brakes on opposite casters $\frac{1}{2}$ With 12 $\frac{1}{2}$ Trendelenburg (TR) & Reverse Trendelenburg (RTR) position Infant Bed, Hospital Baby Crib exporters, Hospital Baby Crib manufacturers