

**Product Name :**  
Semi Fowler Children Bed**Product Code :**  
NMCBS114003**Description :**

Semi Fowler Children Bed is a Semi-Fowler Bed with mechanically operated Back rest by crank mechanism. Size: 1880 x 880 x 580 mm. Epoxy coated mild steel frame work with 2 section top and S.S. side railings with Up & Down movement. With S.S. Shelf at the bottom for keeping blankets, basin etc. Tubular Head & Foot Bows also made of stainless steel. Mounted on four silent swivel castor (all with brakes). Provision for I.V. Rod. Freight saving knock-down construction. Type: With castor / Without castor. Semi Fowler Children Bed, Hospital Baby Crib exporters, Hospital Baby Crib manufacturers