



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Infant Bed Child Cot S.S Adjustable Height

Product Code :
NMCBS114006



Description :

Infant Bed Child Cot S.S Adjustable Height Infant Bed / Cot. Height adjustable approx. 850 - 1000 mm. Completely made of stainless steel. 4 swivel castor (all with brakes). Size 750x400x850/1000mm
1000x400x850/1000mm Infant Bed Child Cot S.S Adjustable Height, Hospital Baby Crib exporters, Hospital Baby Crib manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,