

Product Name :
Anti Bedsore Mattress**Product Code :**
NMCBS115002**Description :**

Anti Bedsore Mattress
Anti Bedsore Mattress based on alternative pressure pad & pump system
Power supply: 220V/50Hz, 1A/250V, Air output: 6-8 liter/min
Pressure range: 30 mmHg - 100 mmHg
Top cover: Nylon / PU with button
Size: 190cm(L) x 85cm(W) x 11cm (H)
Exchangeable air cell, water proof.
Low Noisy, Long life time benefit
Size: Without CPR valve / With CPR valve
Anti Bedsore Mattress, Hospital Beds Mattress exporters, Hospital Beds Mattress manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,