

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Anti Bedsore Mattress

**Product Code :**  
NMCBS115002



**Description :**

Anti Bedsore Mattress is an Anti Bedsore Mattress based on alternative pressure pad & pump system. Power supply: 220V/50Hz, 1A/250V, Air output: 6-8 liter/min. Pressure range: 30 mmHg to 100 mmHg. Top cover: Nylon / PU with button. Size: 190cm(L) x 85cm(W) x 11cm (H). Exchangeable air cell, water proof. Low Noisy, Long life time benefit. Size: Without CPR valve to With CPR valve. Anti Bedsore Mattress, Hospital Beds Mattress exporters, Hospital Beds Mattress manufacturers

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)