

Product Name :
Height Adjustable Over Bed Table**Product Code :**
NMCBS118002**Description :**

Height Adjustable Over Bed Table We are engaged in providing our clients with Height Adjustable Over Bed Table that are available in different models and sizes at market leading prices. These products are widely appreciated by the customers for their corrosion & abrasion resistance, excellent performance, low maintenance and longer functional life. Furthermore, our products can be customized as per the specifications laid down by the clients. Other information : Over bed table with Laminated Wooden Top Size: 750 \times 350 \times 820 / 1050 mm Framework made of Stainless steel Adjustable height from 820 - 1050 mm by gas spring 4 swivel castor (2 with brake) Height Adjustable Over Bed Table, Hospital Bedside Tables exporters, Hospital Bedside Tables manufacturers