

**Product Name :**  
Electric Multifunction Obstetric Bed**Product Code :**  
NMCBS120004**Description :**

Electric Multifunction Obstetric Bed This table is designed for women's obstetric and gynaecological operation, Easy operation, safe and reliable, economical and quick ½ Operating positions - up / down, adjustable back rest & reverse trendelenburg position ½ It adopts linear bearing to control the telescopic working table ½ Size: 2000/1450 ½ 900/660 ½ 620/920 mm ½ All motion controlled with T-Motion Power Device (Actuator) & Battery backup ½ Adjustable leg support and sliding grip holder ½ Main framework made of Epoxy Coated steel ½ High quality beautiful ABS engineering plastic head panel ½ Two easy lifting guardrails (1 on each side), which are safe and reliable, and can be fixed upward and downward ½ 125mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press ½ Equipped with 8cm thick high quality foam mattress Electric Multifunction Obstetric Bed, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers