

Product Name :
Delivery Bed**Product Code :**
NMCBS120007**Description :**

Delivery Bed Size: 1800L x 550W x 800H mm
Adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips
Supports for leg holders with solid mechanism to fix leg holders at required height and width
Supplied with 2 knee crutches, adjustable in height and width
Framework made of mild steel epoxy coated finish with upholstered top
Supplied with S.S. bowl
Knockdown construction Delivery Bed, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers